

TITLE _____

Form Page No. _____

Post dialysis volumes & observation

		<u>Vol</u>	<u>Conc'd to</u>
65-75	Quite turbid	38 ml	6.5 ml
42-50	Sl. turbid	30 ml	~ 4.5 ml
Buffer	Clear	36 ml	~ 4 ml

Concentrate to ~ 5 ml using YM 10 membranes

Then spin turbid ones & check sups vs pellets on gels - to decide whether to include pellets.

65-75	} Filtered thru 0.45 μ Acrodisc / Millar type filters	→ Filter burst; recovered 1.4-1.5 ml
42-50		→ Recovered ~ 2 ml

Note both were contaminated w/ red-like bacteria

To Helen Hickman (for Abs, Inc) Redacted
~ 1.4 ml 65-75 above; to be sent on dry ice
3 rabbits200 μ g injections at 1, 7, 21, 35, 56 days
Bleeds at 28 & 63 daysEstimate 1 μ g/ml (see next page)

... enough for first 2 injections

will prepare more for subsequent ones

To Page No. _____

Witnessed & Understood by me, _____

Date _____

Investigated by _____

Date _____

Recorded by _____

Redacted