

Form Pag No. _____

Redacted

Preparation of material ^{from} S-200 in Sarkosyl for ad production

To be sent to AB, Inc.

Helen Hickman:

3 rabbits

→ 1, 14, 28, 35 days? or all at once?

Initially 5 injections ea after adjuvanting (mix w/
2 blocks, at 31 & 61 days equal vol)

Can be 1 ml per injection (15 ml total)

1 mg per rabbit (total)

Should give her 0.5 mg extra, for titrating

Note that this will titrate abs to contaminants also.
(plus binding)

Likely that Sarkosyl will interfere with adjuvanting.

Will make pool, & dialyze vs 20 mM Tris-HCl, pH 8

Then concentrate

Do likewise for

Pool from another part of column } as controls
BufferPools: 65-75 C.I. pool; excluding contaminants as
much as possible

Vol = ~31 ml

Control 42-50 ~26 ml

Buffer: 30 ml

Dialyze each vs 1 L 3X 50 mM 1/9 → 9 AM 1/10
RT buffer at 15° → 4°C → 9 AM 1/11

→ 9 AM 1/13 (vs 500 ml)

To Pag No. _____

Witnessed & Understood by me,

Date

Inv nt d by

Dat

Rec rd d by

KEC

Redacted