

Form Pag No. \_\_\_\_\_

Redacted

Preparation of material <sup>from</sup> S-200 in Sarkosyl for ad production

To be sent to AB, Inc.

Helen Hickman:

3 rabbits

→ 1, 14, 28, 35 days? or all at once?

Initially 5 injections ea after adjuvanting (mix w/  
2 blocks, at 31 & 61 days equal vol)

Can be 1 ml per injection (15 ml total)

1 mg per rabbit (total)

Should give her 0.5 mg extra, for titrating

Note that this will <sup>(high</sup> titer abs to contaminants also.  
binding)

Likely that Sarkosyl will interfere with adjuvanting.

Will make pool, &amp; dialyze vs 20 mM Tris-HCl, pH 8

Then concentrate

Do likewise for

Pool from another part of column } as controls  
BufferPools: 65-75 C.I. pool; excluding contaminants as  
much as possible

Vol = ~31 ml

Control 42-50 ~26 ml

Buffer : 30 ml

Dialyze each vs 1 L 3X 50 mM 1/9 → 9 AM 1/10  
RT buffer at 15° → 4°C → 9 AM 1/11

→ 9 AM 1/13 (vs 500 ml)

To Pag No. \_\_\_\_\_

Witnessed &amp; Understood by me,

Date

Inv nt d by

Dat

Rec rd d by

KEC

Redacted